

PATIENT INSTRUCTIONS FOR RIGISCAN

The device will be supplied with new batteries. There will be 1 day worth of additional batteries inside of the carrying case.

Do not turn on the Rigiscan until loops have been placed on the penis.

FIRST NIGHT AT BEDTIME:

- Place base loop around the base of the penis.
- Place tip loop around the tip of the penis.
- If loops are too small, insert fingers in loops and slowly widen. (The device must be turned off when doing this).
- If the loops are too large, the Rigiscan will adjust them to size when the unit is turned on.
- Turn the device on by pushing the on/off button.

UPON AWAKENING:

- Turn off the device: remove the loops.
- Change the batteries; discard or recycle used batteries immediately.
- Place the new batteries in the Rigiscan and store the device in the carrying case.

SECOND NIGHT AT BEDTIME:

- Repeat above process

Return the Rigiscan to Michigan Institute of Urology office the morning after the second exam.

REMEMBER:

- Do not turn the device on until the loops have been placed on the penis as described above.
- Do not take sleeping medications or muscle relaxants.
- Take all daily medication as prescribed (except sleeping medication and muscle relaxants).
- Do not use alcohol or other drugs or drink coffee, tea or colas after 3 pm.
- Do not attempt sexual intercourse (including oral sex or masturbation) while using the Rigiscan.
- Try to empty your bladder before using the device: if you need to urinate, turn the device off without removing the loops and then urinate. You may leave the device off for up to 10 minutes. ***Remember to turn the device back on before returning to bed.***
- Always place new batteries in the device for the second night of use. One set will be in the device, and a second set will be in the carrying case.

THE RETURN DATES ARE AS FOLLOWS. PLEASE HAVE THE EQUIPMENT RETURNED BY 9 AM ON THE DUE DATE:

If you receive the Rigiscan:

Monday—return by 9 am Wednesday

Thursday—return by 9 am Monday

Tuesday—return by 9 am Thursday

Friday—return by 9 am Monday

Wednesday—return by 9 am Friday